


C-SMART Pilot

True wireless Endo motor
with Apex locator



Appropriate slope to
ensure clear observation
from all angles



5 inch **Colorful** Screen




7 mainstream bands

30+ Endo file systems

200+ Endo file parameter

Multi-Function Buttons





Imported Motor
Powerful · Mute · Longlife

Speed

**150-1000
RPM**

Torque


**0.6-3.9
N.CM**

Apex locator built in the handle


Just plug in test wire and use,
get rid of the main machine,
get rid of inconvenient



Integrated insulation
Contra angle



Built in root
measurement
buckle



Old one

Built in LED

300° Rotation angle

